

PIONEERING ADVANTAGES



- Made in Germany with patented technology
- Ideal solution for medical and aesthetics
- Applied whole body on skin and mucosal
- Safe, Fast, Easy and Effective
- Varieties of 24K gold-plated applicators
- Medical grade pure N2O liquid work at -89°C
- Sealed system only negligible leaking during storage

SM+SC

翔林國際有限公司

Shianglin Medical System Co.

100台北市中正區和平西路2段141號3樓之3

Tel 02-23123179 Fax 02-23122711



CRYOALFA FROM GERMANY


CRYO-SURGICAL PEN

IDEAL SOLUTION FOR MEDICAL AND AESTHETICS

LUX SERIES




Lux Standard + Lever Control



Lux Derma



Lux Contact



Lux Gyn

LUX MEGA BOX


- Lux Standard
- Lux Lever Control
- Lux Derma 4mm
- Lux Contact 5mm
- 16g Cartridge
- 25g Cartridge
- Gift Box



SUPER SERIES



Super Standard



Super Contact

CARTRIDGE



16g N₂O cartridge
(approx. 50 exams)



25g N₂O cartridge
(approx. 78 exams)

LUX GYN BOX

- Lux Standard
- Lux Lever Control
- Lux Gyn 9cm
- 25g Cartridge
- Gift Box




TREATMENT




Skin Cancer



Keloid



Wart



Mole

Other Applications:

Hoemangioma, Fibromas, Polyp, Skin tag, Aging pigmentation, Lentigines, Acne nodulus, Juvenile Warts, Keratosis actinica, Keratoma senile, Basilioma, Naevus pigmentosus, Ulcer . . .

CLINICAL DEMONSTRATIONS

YouTube

CRYOALFA

For details please contact us:

2544 2655 / 9495 4478 / 9626 2544

<http://bit.ly/cryoalfa-demo>

