

PORT-X

PORTABLE DENTAL X-RAY

Cat. No
X4000

The PORT-X IV provides total portability combined with the latest digital technology. Ideal for both dental and extremity radiographs with vet specific icons and pre-programmed smart exposure settings.

The Port-X IV combines the latest portable technology with the highest specifications amongst portable dental x-ray generators combined with built in safety features including a lead lined tank assembly and tube.

- Light and small, enabling total portability and use in any location
- Small dose of radiation, ensuring patient and operator safety
- High performance tube and HFG, providing sharp and clear images
- Designed for use with both digital sensors and CR image plates
- Supplied with remote operating switch and scatter shield
- Fixed kV (70) and mA (2) means the only variable is the exposure time
- X-ray tube manufactured by Toshiba
- PORT-X IV can scan up to 400 times * with one full charge of the built-in battery. The SLEEP MODE function adjusts the amount of power consumed by the device for efficient battery management.
- back scatter shield (optional extra)

* Exposure time: 0.5 second

National radiology guidelines vary on the use of hand-held generators. Please contact us for availability in your country.

PORT-X SUPPLIED WITH:

- » Tripod
- » Remote operating switch

1.5 kg



SPECIFICATIONS

WEIGHT	DIMENSIONS	TUBE VOLTAGE	TUBE CURRENT
1.5k g	140 x 173 x 254 mm	70 kV	2 mAs
FOCAL SPOT	EXPOSURE TIME	FDA & CE APPROVED	BATTERY
0.4mm	0.05 – 1.6 sec	✓	Li-Polymer 22.2V DC 1,000mAh/25C