

NEW

petMAP *graphic II**



The next generation petMAP *graphic* for all your BP/HR measurement needs

petMAP *graphic II* provides:

- Proven petMAP technology with PPO (PetMAP Proprietary Optimizations). Species and cuff site settings to improve accuracy.
- Easy to use “no menu” touch screen display
- Real-time display of oscillometric “envelope”
- Selectable BP cycle times for measurement reminders
- Tabular display of readings
- Graphic display of last 45 readings and NSV
- SD memory card to record all BP/HR values and “power off” image
- Optional pole mount, cover and AC/DC adapter

Portability and Accuracy like no other!

CONTACT US TODAY TO PLACE YOUR ORDER

(800) 231-6370

(813) 289-5555

Take advantage of our 30 day satisfaction guarantee.

For more information visit www.petmap.com

Manufactured,
Distributed and
Serviced by:



www.cardiocommand.com



**RAMSEY
MEDICAL INC**

Developed by: **RAMSEY
MEDICAL INC**
* Trademark of Ramsey Medical, Inc.

petMAP graphic II

petMAP *graphic II* (g II) is a small, lightweight and portable blood pressure measurement device designed for veterinary use. It is capable of measuring blood pressure (systolic, diastolic and mean arterial pressures) and heart rate. The measured values are displayed and trended graphically and numerically on the device. Data trended on the petMAP *g II* can optionally be saved on an SD memory card via the device's integrated SD memory card slot.

A user interacts with petMAP *g II* using its touch screen, just like many smart phones. The user can set alarms for systolic blood pressure and can set a BP cycle time for automatic reminders to start another BP measurement.



7300 petMAP *graphic II*

Blood Pressure Range: Systolic BP: 40 - 265 mmHg
MAP: 30 - 200 mmHg
Diastolic BP: 25 - 195 mmHg

Heart Rate Range: 40 - 240 bpm

Dimensions: 6.0"H x 3.2"W x 1.25"D

Weight: 1 lb (with batteries)

Battery: Four AA Alkaline (installed)

Available Accessories

Pole Mount
Item #: 9011



Cover
Item #: 9015



AC/DC Adapter
Item #: 9012



International
Adaptor Plugs
Item #: 9008



United Kingdom

Europe

Australia

CritterCuffs™



CritterCuffs included with petMAP *gII* System



*Sold separately, not included in petMAP *gII* System

Blood Pressure Cuffs
Seven sizes included with petMAP *gII*

Size	Item no.
2.0cm	8050
2.5cm	8051
3.0cm	8055
3.5cm	8052
4.0cm	8056
4.5cm	8053
5.5cm	8054
5.0cm	8067*
6.5cm	8057*
8.0cm	8058*
10.0cm	8059*
13.0cm	8060*

RamseyCuffs™



*Sold separately, not included in petMAP *gII* System